

Face Shield Specification (PPE)

Doc No.:		FN/PTSP-FS01010106-FN	Revision &Date:	Rev.01 2020-10-09
Brand		FINESS		
Product Name		Face Shield		
Product Code		FN-FS01010106		
Specification		33x22cm		
Color		Transparent visor+blue		
Material	Visor	PET		
	Headband	Urethane elastic fibre		
	Sponge	Polyurethane foam		
	Button	Resin		
		<p>Product Drawing</p> <p>Front side</p> <p>Back side</p> <p>Sponge: LxWxT 24x3x3cm</p> <p>1.5cm 1.5cm 3.2 x 1.5cm Headband 4x3.3cm Label 18.3cm Visor</p>		
Detailed Item Properties			Criteria	Ref.Standard
Size(cm)	Visor:LxW	(33 x 18)±1	N/A	N/A
	Label:LxW	(33 x 4)±1	N/A	N/A
	Headband:LxW	(32±1) x (1.5±0.5)	Width:≥10mm	EN166:2001
	Sponge:LxWxT	(24±1)x[(3x3)±0.5]	N/A	N/A
Thickness	Visor:25 si		N/A	
Transmittance	≥74.4%		≥74.4 %	EN166:2001
Diffusion of light	≤0.5(cd/m ² · lx)		≤0.5(cd/m ² · lx)	EN166:2001
Stability	Stable at an elevated temperature		no apparent deformation	EN166:2001
Protection	Protection against liquid splashes		EN166:2001	
Fogging Resistance	Free from fogging		EN166:2001	



FINESS HEALTHCARE GROUP CO.,LTD

Product Technical Specification

Intended Use	Used to protect the wearers from infectious fluids or harmful biological agents splash or spray onto face.
Classification	Category III according to REGULATION (EU) 2016/425 Annex I
Remarks	Non-sterile, Single-use only. The front and back side of face shield are covered with protective films, the film labels should be removed before use.
Packaging Information	Carton Size: 64x35x38cm, 1pc/bag, 200pcs/CTN Gross Weight: 7kg The carton size and gross weight will be varied in accordance with the actual situation.

Document Revision History:

Ver.	Revision Content Description	Effective Date	Complied by	Reviewed by	Approved by
00	Initial Version	2020-04-20	Felix Duan	Ruby Fu	Michelle Zhang
01	Anti-fogging property updated	2020-10-09	Felix Duan	Ruby Fu	Michelle Zhang